



QView™ Head positioning Device (HPD)

E.T. tube support and viewing mirror for use with prone patients



- Gain clear view and access to patient's mouth, nose, eyes and tracheal tube through the slotted mirror design.
- Comfort and support patient's head and neck in the contoured cradle.
- Helps control infections and eliminate reprocessing with single-use.

The Q View Head Positioning Device (HPD) is for use with patients who need to be placed in a PRONE position for an operation. It is designed to support not only the head, but also ensure that ET Tube is positioned and supported and also the patient's face can be monitored by the clinician via the viewing mirror. The Q View is lightweight and latex free.



Extra care for sensitive skin:

A layer of 3M biocompatible material is added onto the skin contacting surface, giving extra care to the patients with sensitive skin.

Product Number:	QPHD-Adult
Product type:	Single use
Intended use:	Adult
Color :	Natural
Dimension:	H x W x L: 15.0 x 24.5 x 30.5 cm
Cushion Material:	Polyurethane Foam
Skin Contacting Surface:	3M Nonwoven Medical Tape #1776
Mirror size:	H x W x L: 0.165 x 20.0 x 25.0 cm
Mirror material:	Acrylic
Packaged:	Non sterile, 6 pcs per box
Storage:	0-30 °C, in box. Keep away from direct sun light

Shenzhen Envisen Industry Co. Limited

2nd Floor, Block 1, 40 Jianlong Street, Baoan Community, Heng Gang Town,
Long Gang District, Shenzhen, China 518115

Ph.: + 86 755 88827227; Fax: +86 755 28604953; Email: info@envisen.com; Website: www.envisen.com



Qview Specification List

Product Code	QView	QView - 3	QView - P	QView - 3S	QView - PS
Dimension (L x W x H) cm	30.5 X 24.5 X 15	30.5 X 24.5 X 12	30.5 X 24.5 X 12	30.5 X 24.5 X 12	30.5 X 24.5 X 12
Cushion Material	Polyurethane Foam	Polyurethane Foam	Polyurethane Foam	Polyurethane Memory Foam	Polyurethane Memory Foam
Skin Contact Surface	3M Micro Foam	3M Micro Foam	3M Micro Foam	N.A.	N.A.
Mirror Size (L x W x H) cm	25 X 20 X 0.2				
Packaging	6 pcs per box	8 pcs per box	8 pcs per box	8 pcs per box	8 pcs per box
Storage	0-30°C in box, avoid light exposure				
Supplementary Accessory	None	None	Compatible with ProneView® Protective Helmet	Supported on QRack	Compatible with ProneView® Protective Helmet