



To select a correct BRIDGE™ SpO<sub>2</sub> Sensor reference number...

Pick up a sensor type (Prefix) from below,



CF

Finger Clip

(Patient Weight >20kg)



CFS

Small Finger Clip

(Patient Weight >20kg)



URL

Ultra Soft Finger

(Patient Weight >20kg)



URS

Small Ultra Soft Finger

(Patient Weight >20kg)



EC

Ear Clip w. hanger

(Patient Weight >30kg)



Y

Multi Site Sensor

(Patient Weight > 1kg)



W

Wrap Sensor

(Patient Weight >1kg)



Prefix

Suffix


			-											
--	--	--	---	--	--	--	--	--	--	--	--	--	--	--



SpO<sub>2</sub> sensors' Suffix and Cables ordering information please refer to the following pages




















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Agilent Technologies (See Philips Medical)									
Armstrong Medical									
AD-900, 1000, 2000, 3000, 4000	BCI 3044	<b>1203-120</b>	120	C03  	N/A	X-126033-240	240	C60  	C33  
Artema									
MM200	Not Available	<b>6012-295</b>	295	C12  	N/A	Direct Connect	N/A	N/A	N/A
Diascope	Nellcor DS-100A	<b>1105-90</b>	90	C05  	59831-0	X-606016-240	N/A	C60  	C16  
S&W	S&W SL033057	<b>6015-295</b>	295	C15  	N/A	Direct Connect	N/A	N/A	N/A
Baxter									
ASAT	SIMED 600015	<b>6123-295</b>	295	C23  	N/A	Direct Connect	N/A	N/A	N/A
Explorer/Monarch	Nellcor DS-100A	<b>1105-90</b>	90	C05  	EC-8	X-116005-240	240	C60  	C05  
BCI (Biochem Int'l)									
3100, 3101, 6100, 9100, Advisor, Autocorr, MiniTorr, Mini Torr Plus	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI 3109	X-126033-240	240	C60  	C33  
Same as Above	Not Available	<b>1233-295</b>	295	C33  	N/A	Direct Connect	N/A	N/A	N/A
3300, 3301, 3302, 3303, 3304, 3401, 3402	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI3311L	X-126003-240	240	C60  	C03  
Same as Above	Not Available	<b>1203-295</b>	295	C03  	N/A	N/A	N/A	N/A	N/A























## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Biolight (Guangdong) Meditech Co., Ltd.										
BLT series M69, M66, M6, M12	Not Available	<b>4280-295</b>	295	C80 	N/A	N/A	N/A	N/A	N/A	
Biosys Co., LTd.										
BPM-200, Vital Sign Monitor SENTRY	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 	
Bird										
4400	CSI 511-10L	<b>1831-295</b>	295	C31 	N/A	Direct Connect	N/A	N/A	N/A	
Same as above	CSI 511SD	<b>1808-90</b>	90	C08 	CSI 518LD	X-186031-240	240	C60 	C31 	
					CSI 940HP	X-186008-240	240	C60 	C08 	
Bruker (See Odam)										
Burdick (See Spacelabs Burdick)										
CAS Medical Systems Inc.										
OscilloMate 9002, 9303, 9302S	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8500VI	X-136004-240	240	C60 	C04 	
OscilloMate 9002, 9303, 9302S	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8500I	X-136004-120	120	C60 	C04 	
Choice Electronic Technologie Co.										
MMED6000DP, F3, F4, F5, F6, F7, md300-D	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03 	N/A	Direct Connect	N/A	N/A	N/A	
MMED6000DP-SF, M3, M6	Not Available	<b>4080-295</b>	295	C80 	N/A	Direct Connect	N/A	N/A	N/A	

























## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Colin									
Press-Mate SAT, BP88, BP306, BP508, BP808, BP88S	Nellcor DS-100A	1105-90	90	C05 	SCP-10/MC-10	X-116041-295	300	C60 	C41 
					EC-8	X-116005-240	240	C60 	C05 
Corometrics (See GE Healthcare - Corometrics) CSI (Criticare System Int'l)									
501, 503, 504, 507, 508, 581, 1100, 8100, Poet (602-3A), Scholar	CSI 511SD	1808-90	90	C08 	CSI 518LD	X-186031-240	240	C60 	C31 
					CSI 940HP	X-186008-240	240	C60 	C08 
Same as Above		1808-295	295	C08 	N/A	Direct Connect	N/A	N/A	
Same as Above	CSI 511-10L	1831-295	295	C31 	N/A	Direct Connect	N/A	N/A	
503 Spot, 504 Plus, 504DX, 506DXN	CSI 511SD	1808-90	90	C08 	CSI 940HP	X-186008-240	240	C60 	C08 
					N/A	Direct Connect	N/A	N/A	N/A
Same as Above		1808-295	295	C08 	N/A	Direct Connect	N/A	N/A	
Critikon (See GE Healthcare - Critikon)									
Datascope									
Accusat, Accutorr 3/4 SAT, Accutorr Plus, (Datascope SpO <sub>2</sub> ), Passport (Datascope SpO <sub>2</sub> )	Datascope 0600-00-0026-01	1909-90	90	C09 	0012-00-1099 02	X-196051-295	295	C60 	C51 
Same as Above	Datascope 0020-00-071-01	1951-295	295	C51 	N/A	Direct Connect	N/A	N/A	













































## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Datascope *Continue										
Expert, Passtport (Nellcor SpO <sub>2</sub> )	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 	
Same as Above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A	
Accutorr Plus (Nellcor SpO <sub>2</sub> )	Nellcor DS-100A	<b>1105-90</b>	90	C05 	SCP-10/ MC-10	X-116041-295	300	C60 	C41 	
Digicare										
Digimax 5500	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03 	N/A	X-126034-150	150	C60 	C34 	
					BCI3311L	X-126003-240	240	C60 	C03 	
Dima Italia										
All	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 	
Same as Above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A	
Dolphin Medical										
210, 2100, 2150 or Dolphin ONE Module 701 Technologies employed	Dolphin 201, 360 digital sensors	<b>4170-50</b>	50	C70 	110	X-416041-240	240	C60 	C41 	
					N/A	X-416070-240	240	C60 	C70 	
Same as Above	N/A	<b>4141-295</b>	295	C41 	N/A	Direct Connect	N/A	N/A	N/A	






















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
<b>Dräger</b>									
Draeger Oxipac 2000, 2500	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI3311L	X-126003-240	240	C60  	C03  
Same as Above	Not Available	<b>1203-295</b>	295	C03  	N/A	Direct Connect	N/A	N/A	N/A
Vialert 1000, Cicero, GammaXL, DeltaXL, Kappa	Nellcor DS-100A	<b>1105-90</b>	90	C05  	N/A	X-326039-240	240	C60  	C39  
Narkomed	Nellcor DS-100A	<b>1105-90</b>	90	C05  	Pre-Amp	X-116038-240	240	C60  	C38  
Infinity Delta series, Kappa Series (w. Infinity MultiMed 12 Pod, Nellcor Modules)	Nellcor DS-100A	<b>1105-90</b>	90	C05  	Not Available	X-326073-150	150	C60  	C73  
<b>EDAN Instrucments Inc.</b>									
M3, M3A, M8A, M8B, M9 series (Edan SpO <sub>2</sub> Module)	EDAN SH1	<b>2935-295</b>	295	C35  	N/A	Direct Connect	N/A	N/A	N/A
<b>Elmed</b>									
PSO1000	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI3311L	X-126003-240	240	C60  	C03  
Same as above	Not Available	<b>1203-295</b>	295	C03  	N/A	Direct Connect	N/A	N/A	N/A
<b>Entri Medical</b>									
Propaq 102, 104, 106	Nellcor DS-100A	<b>1105-90</b>	90	C05  	Pre-Amp	X-116038-240	240	C60  	C38  




















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
EnviteC-Wismar GmbH									
PulsiQuant 8040	F-3222-12	<b>1203-120</b>	120	C03 	X-4222	X-126003-240	240	C60 	C03 
PulsiQuant P	F-3212-9	<b>1105-90</b>	90	C05 	X-4224	N/A	240	C60 	C34 
Fukuda Denshi									
HG302 Module, DS-5100E, DS-5300W	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 
Same as above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A
LS-100	Not Available	<b>6244-33</b>	33	C44 	N/A	Direct Connect	N/A	N/A	N/A
LX-5230	Not Available	<b>6245-33</b>	33	C45 	N/A	Direct Connect	N/A	N/A	N/A
Gmabro									
BP806	Nellcor DS-100A	<b>1105-90</b>	90	C05 	N/A	Direct Connect	N/A	N/A	N/A
GE Healthcare - Datex-Ohmeda									
AS/3; CS/3, CH/S; CH-RS, CG/2S; Cardiocap I, II, 5; Capnomac Ultima; Light Monitor; M-Estpr; Oscar II; OSP-200; Satellite Trans Plus, ULTSV I; ULTSVO; ULTSIO; ULTS00	Datex SAS-F	<b>1501-90</b>	90	C01 	OXY-C3	X-156014-295	295	C60 	C14 
Same as Above	OXY-F-DB	<b>1501-200</b>	200	C01 	OXY-C3	X-156014-295	295	C60 	C14 
Same as Above	OXY-F4-N	<b>1514-295</b>	295	C14 	N/A	Direct Connect	N/A	N/A	N/A
















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
GE Healthcare - Datex-Ohmeda *Continue									
AS/3, CS/3, S/5 with M-NSAT Module (Nellcor Non-OxiMAX)	Nellcor DS-100A	<b>1105-90</b>	90	C05 	SCP-10 MC-10	X-116041-295	295	C60 	C41 
i4 with M-NSAT Module	Nellcor DS-100A	<b>2679-90</b>	90	C79 	DOC-10	X-266968-295	295	C69 	C68 
S/5 FM Port	OXY-F4-GE/ OXY-E4-GE	<b>1564-295</b>	295	C64 	N/A	Direct Connect	N/A	N/A	N/A
Ohmeda 3700, 3710, 3770, 3775, Biox 3740, 3760; 3800, 3900, 4700, TuffSat, OxiCap, 5250RGM, Reascal II	OXY-F4-H	<b>1422-295</b>	295	C22 	N/A	Direct Connect	N/A	N/A	N/A
Same as Above	N/A	<b>1401-90</b>	90	C01 	N/A	X-146022-240	240	C60 	C22 
Same as Above	OXY-F1-H	<b>1422-100</b>	100	C22 	N/A	Direct Connect	N/A	N/A	N/A
GE Healthcare - Marquette Hellige									
7000 series, TRAM A.R., TRAM 200SL, 400SL, 600SL, 800A, 800SL, SMU, SMK, SMV	OXY-F4-H	<b>1422-295</b>	295	C22 	N/A	Direct Connect	N/A	N/A	N/A
Same as above	Not Available	<b>1401-90</b>	90	C01 	N/A	X-146022-240	240	C60 	C22 
Eagle, Dash, SOLAR, TRAM Modular, 250SL, 450SL, 650SL, 850A, 850SL	Nellcor DS-100A	<b>1105-90</b>	90	C05 	223 394 01	X-316050-240	240	C60 	C50 
Hellige Eagle 1000, 3000, SMU, SMK, SMV	Nellcor DS-100A	<b>1105-90</b>	90	C05 	303 443 58	X-316055-295	295	C60 	C55 


























## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
GE Healthcare - Marquette Hellige *Continue									
DASH 2500/3000/4000/5000 Solar 8000/8000Mi/9500 (Tram 451N5/851N5) Corometrics 250	Nellcor DS-100A	<b>2679-90</b>	90	C79 	Unknown	Upon Request	295	C69 	Waiting Confirm 
DASH 3000/4000/5000 Solar 8000/8000Mi/9500 (Tram 451M/851M) Corometrics 259 (M Technology Employed)	Unknown	<b>1766-90</b>	90	C66 	Unknown	Upon Request	295	Waiting Confirm Photos Not Available	Waiting Confirm 
GE Healthcare - Corometrics									
Corometrics 188, 129, 129F, 128, 128F	Nellcor DS-100A	<b>1105-90</b>	90	C05 	223 394 01	X-316050-240	240	C60 	C50 
					EC-8	X-116005-240	240	C60 	C05 
Corometrics 120, 510, 511, 118	Nellcor DS-100A	<b>1105-90</b>	90	C05 	N/A	Upon Request	180	C60 	C27 
Corometrics 556, 570	OXY-F4-H	<b>1422-295</b>	295	C22 	N/A	Direct Connect	N/A	N/A	N/A
GE Healthcare - Critikon									
Critikon Dinamap Oxytrak, 8280	Critikon 8997	<b>3102-90</b>	90	C02 	N/A	[None Offered]	N/A	N/A	N/A
Critikon Dinamap 8700/ 8710/ 8720 PLUS	Critikon 9084 + 8731	<b>3932-295</b>	295	C32 	N/A	Direct Connect	N/A	N/A	N/A
Critikon Dinamap 8700/ 8710/ 8720 PLUS, Oxyshuttle	Critikon 9084	<b>3946-90</b>	90	C46 	8731	X-396132-295	295	C61 	C32 




















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
GE Healthcare - Critikon *Continue									
Critikon Dinamap 9700/ 9710/ 9720 PLUS, Compact, MPS, Select	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 
Same as above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A
Dinamap Pro 100-400, Pro 1000	Nellcor DS-100A	<b>1105-90</b>	90	C05 	SCP-10/ MC-10	X-116041-295	295	C60 	C41 
Dinamap Pro V2	Nellcor DS-100A	<b>2679-90</b>	90	C79 	DOC-10	X-266968-295	295	C69 	C68 
Goldway Industrial Inc.									
UT4000A	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03 	N/A	X-696080-150	150	C60 	C80 
Same as Above	Not Available	<b>6980-295</b>	295	C80 	N/A	Direct Connect	N/A	N/A	N/A
Hellige (See GE Healthcare - Marquette Hellige) Healthdyne									
All	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8500VI	X-136004-240	240	C60 	C04 
Same as above	Not Available	<b>1304-295</b>	295	C04 	N/A	Direct Connect	N/A	N/A	N/A
Hewlett Packard (See Philips Medcal) Innomed									
Oxycard, Innocare-T0	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8500VI	X-136004-240	240	C60 	C04 






















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Innomed *Continue										
Same as above	Not Available	<b>1304-295</b>	295	C04 	N/A	Direct Connect	N/A	N/A	N/A	
Invacare										
All	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03 	BCI3311L	X-126003-240	240	C60 	C03 	
Same as above	Not Available	<b>1203-295</b>	295	C03 	N/A	Direct Connect	N/A	N/A	N/A	
Invivo Resarch										
4500, 4500 Plus, Millenia, Omega 1445, Omni-Trak TVS, Sat Scout	Invivo 9383	<b>6332-295</b>	295	C32 	N/A	Direct Connect	N/A	N/A	N/A	
4500 Plus 2	Nellcor DS-100A	<b>1105-90</b>	90	C05 	N/A	X-636032-240	240	C60 	C32 	
Ivy Biomedical										
405A, 405D	Nellcor DS-100A	<b>1105-90</b>	90	C05 	Pre-Amp	X-116038-240	240	C60 	C38 	
405P, 405T	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 	
Same as above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A	
Koike Medical										
Wrist SO2	Unknown	<b>Upon Request</b>	33	C44 	N/A	N/A	N/A	N/A	N/A	















































## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Kontron										
7840, 7845, Minimon 7137, Micromon, Supermon, Kolormon, Trakmon	Kontron 0608010	N/A	295	C11 	N/A	Direct Connect	N/A	N/A	N/A	
Same as above	Nellcor DS-100A	1105-90	90	C05 	N/A	Upon request	240	C60 	C11 	
Larsen & Toubro Limited (see L&T Medical)										
L&T Medical										
Steller P, Star 50, Galaxy, Planet, Planet 50, Cosmic	BCI 3044/BCI3044N	1203-120	120	C03 	N/A	Direct Connect	N/A	N/A	N/A	
Lynn Medical										
All	BCI 3044/BCI3044N	1203-120	120	C03 	BCI3311L	X-126003-240	240	C60 	C03 	
Same as above	Not Available	1203-295	295	C03 	N/A	Direct Connect	N/A	N/A	N/A	
Magna Medical										
MON-MTEK-30, 33, 35, 40, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, MON-MTEK-MOD2/ MOD3	BCI 3044/BCI3044N	1203-120	120	C03 	BCI3311L	X-126003-240	240	C60 	C03 	
Same as above	Not Available	1203-295	295	C03 	N/A	Direct Connect	N/A	N/A	N/A	
Mallinckrodt/ TYCO/ Covidient										
N-100 (w/ High Sat), N-200	Nellcor DS-100A	1105-90	90	C05 	E610-01	X-116038-240	240	C60 	C38 	
Same as above	Nellcor DS-100A + EC8	1105-295	295	C05 	E610-01	X-116038-240	240	C60 	C38 	











































## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Mallinckrodt/ TYCO/ Covidient *Continue									
NPB-290/5, N-390, N-3000, N-6000	Nellcor DS-100A	<b>1105-90</b>	90	C05  	SCP-10/ MC-10	X-116041-295	295	C60  	C41  
Same as above	Nellcor DS-100A + EC8	<b>1105-295</b>	295	C05  	SCP-10/ MC-10	X-116041-295	295	C60  	C41  
NPB-40	Nellcor DS-100A	<b>1105-90</b>	90	C05  	EC-8	X-116005-240	240	C60  	C05  
N-65, N-85	Nellcor DS-100A	<b>2679-90</b>	90	C79  	DEC-8	X-266979-240	240	C69  	C79  
NPB-500, 595, 600, 650	Nellcor DS-100A	<b>2679-90</b>	90	C79  	DOC-10	X-266968-295	295	C69  	C68  
					DEC-8	X-266979-240	240	C69  	C79  
Marquette (See GE Healthcare - Marquette Hellige) Mediana									
YM-2002 series	Nellcor DS-100A	<b>1105-90</b>	90	C05  	DOC-10	X-266968-295	295	C69  	C68  
					39783	X-266979-240	240	C69  	C79  








## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
<b>MEDLAB</b>									
NANOX2, NANOXeco	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI3311L	X-126003-240	240	C60  	C03  
PEARL 10, NANOX 10, NANOX10C	Nellcor DS-100A	<b>1105-90</b>	90	C05  	N/A	Direct Connect	N/A	N/A	N/A
CAPNOX, VITRO, elan, POX 10, P-OX100L	Nellcor DS-100A	<b>1105-90</b>	90	C05  	309	X-116034-240	N/A	C60  	C34  
Same as above	Not Available	<b>N/A</b>	295	C34  	N/A	Direct Connect	N/A	N/A	N/A
<b>Medsource</b>									
BPM-200	Nellcor DS-100A	<b>1105-90</b>	90	C05  	EC-8	X-116005-240	240	C60  	C05  
Same as above	Not Available	<b>1105-295</b>	295	C05  	N/A	Direct Connect	N/A	N/A	N/A
<b>M.D.E.</b>									
Escort 100 OPT11A, Scort 300A	Nellcor DS-100A	<b>1105-90</b>	90	C05  	E610-01	X-116038-240	240	C60  	C38  
Escort II POT30; Prism	Nellcor DS-100A	<b>1105-90</b>	90	C05  	EC-8	X-116005-240	240	C60  	C05  
Same as above	Not Available	<b>1105-295</b>	295	C05  	N/A	Direct Connect	N/A	N/A	N/A
<b>MEK</b>									
MP-100, MP-110, MP-111, MP-200, MP-400, MP-500, MP-1000	Not Available	<b>Upon Request</b>	295	C06  	N/A	Direct Connect	N/A	N/A	N/A





















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Mennen										
Horizon XL (Ohmeda SpO <sub>2</sub> )	Mennen 0800-060-540	<b>1422-295</b>	295	C22 	N/A	Direct Connect	N/A	N/A	N/A	
Same as above	Not Available	<b>1401-90</b>	90	C01 	Not Available	X-146022-240	N/A	C60 	C22 	
Horizon XL (Nellcor SpO <sub>2</sub> ), Mercury 1200, MR 1303	Nellcor DS-100A	<b>1105-90</b>	90	C05 	E610-01/ Preamp	X-116038-240	240	C60 	C38 	
Mindray										
PM-600, PM-6000, PM-7000, PM-8000, PM-9000, MEC-1000, MEC2000, MEC509B (Mindray SpO <sub>2</sub> Module)	Mindray Technologie	<b>1105-90</b>	90	C05 	Not Available	X-286035-240	240	C60 	C35 	
Same as above	Mindray Technologie	<b>2835-295</b>	295	C35 	EC-8	X-116005-240	240	C60 	C05 	
Same as above	Mindray Technologie	<b>2835-295</b>	295	C35 	N/A	Direct Connect	N/A	N/A	N/A	
T5, T8 (Nellcor Module)	Nellcor Technology	<b>2679-90</b>	90	C79 	Unknown	X-286917-295	295	C69 	C17 	
T5, T8 (M Technology)	M series	<b>1766-90</b>	90	C66 	Unknown	Upon Request	295	Waiting Confirm	Photos Not Available	
PM-50	Mindray Technologie	<b>1105-90</b>	90	C05 	N/A	Direct Connect	N/A	N/A	N/A	
Minolta										
Pulsox-3	SR-5C	<b>2145-33</b>	33	C45 	N/A	Direct Connect	N/A	N/A	N/A	




















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Nellcor (See Mallinckrodt) Nihon Kohden									
AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200; OLV-1100A, Lifescope	N.K. TL-101S	<b>2218-295</b>	295	C18 	N/A	Direct Connect	N/A	N/A	N/A
Same as above	Not Available	<b>2201-90</b>	90	C01 	Not Available	X-226018-240	240	C60 	C18 
AL-900P Module; BSM-1102, 9510; OGS-2000	Nellcor DS-100A	<b>1105-90</b>	90	C05 	EC-8	X-116005-240	240	C60 	C05 
Same as above	Not Available	<b>1105-295</b>	295	C05 	N/A	Direct Connect	N/A	N/A	N/A
BSM-4104A, 4114A	Nellcor DS-100A	<b>1105-90</b>	90	C05 	SCP-10/ MC-10	X-116041-295	295	C60 	C41 
BSM-2301K, BSM-2303K, BSM-2351, BSM-2353, BSM-4103, BSM-4113, BSM-5105, BSM-5135, BSM-9501J/K, BSM-9800, OPV-1500, ZS-930P, ZS-940P	N.K. TL-201T	<b>2210-160</b>	160	C10 	JL-900P	Upon Request	295	Upon Request	Upon Request
Nonin									
8500, PamSAT 2500, Avant 9600, Avant 9700, Avant 2120, 9840	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8500VI	X-136004-240	240	C60 	C04 
Same as above	Not Available	<b>1304-295</b>	295	C04 	N/A	Direct Connect	N/A	N/A	N/A
8600, 8700, 8800	Nonin 8000AA-1	<b>1304-90</b>	90	C04 	8604X-10	[None Offered]	N/A	C60 	[None Offered]
					8500I	X-136004-120	120	C60 	C04 















































## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Nonin *Continue									
Same as above	N/A	<b>1304-295</b>	295	C04 	8604X-10	[None Offered]	N/A	C60 	[None Offered]
Same as above (Right Angle Connector)	N/A	<b>1349-295</b>	295	C49 	N/A	Direct Connect	N/A	N/A	N/A
Same as above (Right Angle Connector)	N/A	<b>1349-90</b>	90	C49 	8604X-10	Not Available	N/A	C60 	[None Offered]
					8500VI	X-136004-240	240	C60 	C04 
Novamatrix									
510, 511, 505, 515B/C, 520, 7000, 7100, Oxypleth	Novamatrix 8776-00	<b>1621-295</b>	295	C21 	N/A	X-162621-295	295	C26 	C21 
Same as above	Not Available	<b>1601-90</b>	90	C01 	Not Available	X-166021-240	240	C60 	C21 
500	Novamatrix 8660-00	<b>3621-295</b>	295	C21 	N/A	X-162621-295	295	C26 	C21 
Handheld 512, 513, 710, 715, MARSpO <sub>2</sub> , Model 2001	Novamatrix 9769-00	<b>1607-90</b>	90	C07 	N/A	Direct Connect	N/A	N/A	N/A
Odam									
All	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03 	Not Available	X-126030-240	240	C60 	C30 
Ohmeda (See GE Healthcare - Datex-Ohmeda)									



































## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
<b>PaceTech</b>									
Vitalmax 800 Series (Black face, Nellcor SpO <sub>2</sub> )	PaceTech 4510	<b>1105-90</b>	90	C05  	Not Available	X-666047-240	240	C60  	C47  
Model 520	PaceTech 4510 + 4530	<b>6647-295</b>	295	C47  	N/A	Direct Connect	N/A	N/A	N/A
<b>Palco</b>									
All	Not Available	<b>2543-90</b>	90	C43  	N/A	Direct Connect	N/A	N/A	N/A
All	Not Available	<b>2543-240</b>	240	C43  	N/A	Direct Connect	N/A	N/A	N/A
<b>Philips Medical</b>									
78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B, HP CMS, V24/V26, Viridia Series 50, Code Master XL	Nellcor DS-100A + Agilent M1900B	<b>2417-295</b>	295	C17  	N/A	Direct Connect	N/A	N/A	N/A
Same as above, Viridia A1, A3	Nellcor DS-100A	<b>1105-90</b>	90	C05  	M1900B	X-246017-295	295	C60  	C17  
(Same as above)	M1191A	<b>2419-200</b>	200	C19  	M1940A	X-246217-200	200	C62  	C17  
MP5, MP20, MP30, MP40, MP50, MP60, MP70, MP90, M4529A, Telemetry: M2600A (M2601A), C3, M3/M4, IntelliVue	M1191A	<b>2419-200</b>	200	C19  	M1941A	X-246219-200	200	C62  	C19  
Same as above	Nellcor DS-100A	<b>1105-90</b>	90	C05  	M1943A	X-246019-200	200	C60  	C19  
					M1941A	X-246219-200	200	C62  	C19  






















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Physio-Control									
Lifepak 12	Nellcor DS-100A	<b>1105-90</b>	90	C05  	N/A	Direct Connect	N/A	N/A	N/A
Protocol Systems (See Welch Allyn Protocol)									
RGB Medical Devices									
Omicrom FT, Omicron PREMIUM	BCI 3044/BCI3044N	<b>1203-120</b>	120	C03  	BCI3109	X-126033-240	240	C60  	C33  
Omicrom FT, Omicron PREMIUM	Not Available	<b>1233-295</b>	295	C33  	N/A	Direct Connect	N/A	N/A	N/A
S&W (See Artema)									
Schiller									
CM8	Nellcor DS-100A	<b>1105-90</b>	90	C05  	Not Available	X-686052-240	240	C60  	C52  
Argus TM-7	Nellcor DS-100A	<b>1105-90</b>	90	C05  	Not Available	Upon reuest	240	C60  	C54  
Same as above	Not Available	<b>N/A</b>	295	C54  	N/A	Direct Connect	N/A	N/A	N/A
SensorMedics									
Oxyshuttle I, II	SensorMedics 766721	<b>3946-90</b>	90	C46  	9025	X-396142-240	240	C61  	C42  
Same as above	Not Available	<b>3942-295</b>	295	C42  	N/A	Direct Connect	N/A	N/A	N/A























## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Siemens Medical									
Sirecust 700, 630; Sirem SpO <sub>2</sub> Module	Nellcor DS-100A	1105-90	90	C05 	4533840EH50U	X-306020-240	240	C60 	C20 
SC6002XL, SC7000, SC8000, SC9000XL	Nellcor DS-100A	1105-90	90	C05 	3368433	X-306073-150	150	C60 	C73 
Same as above	Nellcor DS-100A	1105-90	90	C05 	3368433	X-306073-220	220	C60 	C73 
MultiPod Module	Nellcor DS-100A	1105-90	90	C05 	3368433	X-306073-150	150	C60 	C73 
Same as above	Nellcor DS-100A	1105-90	90	C05 	3368433	X-306073-220	220	C60 	C73 
SIMED									
S-50, S100, S-100e	SIMED 600015	N/A	295	C23 	N/A	Direct Connect	N/A	N/A	N/A
Smiths Medical PM Inc. See BCI(Biochem Int)									
Spacelabs Burdick									
PRO2	BCI 3044/ BCI3044N	1203-120	120	C03 	BCI3311L	X-126003-240	240	C60 	C03 

















## Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification					
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)	
Spacelabs Medical										
90419, 90501, 90651A-03	Spacelabs 369083-001	<b>3621-295</b>	295	C21 	N/A	Direct Connect	N/A	N/A	N/A	
90419-01, 90651A-07	Spacelabs 015-0130-00	<b>1621-295</b>	295	C21 	N/A	Direct Connect	N/A	N/A	N/A	
Same as above	Not Available	<b>1601-90</b>	90	C01 	Not Available	X-166021-240	240	C60 	C21 	
90351-D/G, 90465, 90466, 90467, 90489, 90651-07/08, 90651A	Spacelabs 015-0130-00	<b>1621-295</b>	295	C21 	Not Available	X-162625-60	60	C26 	C25 	
Same as above	Not Available	<b>3625-295</b>	295	C25 	N/A	Direct Connect	N/A	N/A	N/A	
90351-0/6, 90465, 90466, 90467, 90651A-08, PC-Express, 90489, IM77	Nellcor DS-100A	<b>1105-90</b>	90	C05 	700-0002-00	X-366024-295	295	C60 	C24 	
90651A-08, 90489, 90466, 90467, 90465, 90351-D/G, 90308, PC Express	Spacelabs 015-0130-00 + 175-0646-00	<b>3624-295</b>	295	C24 	N/A	Direct Connect	N/A	N/A	N/A	
90496 (Nellcor SpO2)	Nellcor DS-100A	<b>1105-90</b>	90	C05 	700.0030.00	X-366048-295	295	C60 	C48 	
90496 (Novamatrix SpO2)	Spacelabs 015-0130-00	<b>1621-295</b>	295	C21 	700.0029.00	X-362648-90	90	C26 	C48 	
Takaoka										
All	Not Available	<b>N/A</b>	295	C36 	N/A	Direct Connect	N/A	N/A	N/A	



# Suffix List

Manufacturer	BRIDGE™ SpO <sub>2</sub> Sensor Specification				BRIDGE™ Cable Specification				
	OEM Models	REF Suffix	Cable Length (in cm)	Connector Number & Photo	OEM Models	REF Number	Cable Length (in cm)	Connector Number & Photo (Patient End)	Connector Number & Photo (Monitor End)
Welch Allyn Protocol									
Propaq 206-EL, 204-EL, 202-EL, 106-EL	Nellcor DS-100A	<b>1105-90</b>	90	C05  	EC-8	X-116005-240	240	C60  	C05  
Weinmann									
Oxycount Mini	BCI 3044/ BCI3044N	<b>1203-120</b>	120	C03  	BCI3311L	X-126003-240	240	C60  	C03  
Same as above	Not Available	<b>1203-295</b>	295	C03  	N/A	Direct Connect	N/A	N/A	N/A