







## Precise Signal Delivery




BRIDGE offers a wide range of reusable and disposable temperature probes for all types of patients and compatible with leading international manufacturers

- Measuring body temperature through body orifice (e.g. esophagus, rectum, ear or urinary bladder) or body surface
- Corresponding sensor type for Adult, pediatric, infant and neonate patients
- Proven accuracy:  $\pm 0.1$  °C (25°C – 45°C)
- Standard phone, HP and Siemens connectors available. Other connectors can be customized
- Compatible to YSI 400 series
- Santoprene™ or PVC medical cables
- ISO 9001, 10993 and CE approved




#### Reusable Rectal/ Esophageal

- Autoclavable (134°C, 100 times)
- Available in 4mmØ, 3mmØ and 2mmØ
- High sensitivity and accuracy
- Robust design




#### Reusable Skin

- Autoclavable (134°C, 100 times)
- Made of non allergic materials
- Available in complete patient sizes
- High sensitivity and accuracy
- Robust design




#### Disposable Rectal/ Esophageal

- PVC cable
- Available in 4mmØ, 3mmØ and 2mmØ
- High sensitivity and accuracy
- 'Clinical Clean' or ETO sterile packaging




#### Disposable Skin

- Innovative design with metal disk tip for better measuring
- Maximizing patient comfort
- Available in 3 sizes
- 'Clinical Clean' or ETO sterile packaging




#### Disposable Tympanic

- Soft rubber plug with cotton tip
- Maximizing patient comfort
- PVC cable
- Available in 2 sizes
- Hygienic packaging




#### Foley catheter with temperature probe

- Silicon tube
- Available in 9 sizes (Fr. 8 ~ 24)
- Monitoring bladder temperature
- ETO sterile packaging




#### Disposable twisted wire sensor

- For making varied medical products
- Twisted PVC cable
- Cable length can be customized
- Hygienic packaging




#### Disposable flexible wire sensor

- For making varied medical products
- Flexible PVC cable
- Cable length can be customized
- Hygienic packaging



#### Disposable Copper wire sensor

- For making varied medical products
- Flexible copper cable
- Cable length can be customized
- Tip size < 1mm
- Tolerance of high voltage (5 kV)



#### Extension and Interface Cables

- Available in different lengths 3-5 M
- Molex connector for disposable sensors
- All cable made with durable material
- Service for customized cable